

Highest confidence in sequence accuracy

Platinum SuperFi II DNA Polymerase



Invitrogen™ Platinum™ SuperFi™ II DNA Polymerase is a hot-start, engineered proofreading DNA polymerase. With >300x fidelity than *Taq* polymerase. Its buffer is specially formulated for primer annealing at 60°C

Highlights

- **Excellent accuracy**—>300x the fidelity of *Taq* polymerase
- **Simplified workflow**—no need for a primer melting temperature (T_m) calculator
- **Increased PCR success**—robust amplification of GC-rich targets, DNA of suboptimal purity, and long sequences
- **Improved automation**—high specificity and benchtop stability for 24 hours after reaction setup

>300x *Taq* DNA polymerase fidelity: for preserving DNA sequence accuracy




Figure 1. Fidelity comparison across commercially available DNA polymerases relative to *Taq* polymerase. 3.9 kb PCR amplicons obtained using different DNA polymerases were fragmented with a MuA transposase. Unique molecular identifiers (UMIs) of 12 random nucleotides were introduced during fragmentation to tag each product individually. After next-generation sequencing, reads from the same UMI family were aligned to call errors. Errors were identified only when present in all reads in the UMI family; otherwise they were discarded as sequencing errors. The fidelity values were normalized to *Taq* polymerase's fidelity.



Applications

- Cloning, subcloning, and site-directed mutagenesis
- Long-range PCR (20–40 kb)
- Template generation for sequencing
- Amplification of DNA of suboptimal purity
- High-throughput PCR
- GC-rich PCR



Available application notes

- Site-directed mutagenesis
- Direct PCR with blood
- Colony PCR
- Multiplex PCR

Download the application notes at
thermofisher.com/platinumsuperfi

Robustness, versatility, and simplicity

Universal annealing temperature at 60°C: calculation of T_m for the annealing step is not required




Figure 2. PCR amplification with high specificity and yield with a universal annealing temperature at 60°C. Primer sets with varying annealing temperatures were used to amplify 12 targets from human genomic DNA (gDNA) with a 60°C annealing temperature.

Universal PCR protocol: co-cycling of multiple assays




Figure 3. Time savings and assay co-cycling enabled by universal PCR protocol. Multiple PCR assays can be cycled together using one protocol with a primer-annealing annealing temperature and the extension time selected for the longest amplicon.

Ordering information

Product	Quantity	Cat. No.
Platinum SuperFi II DNA Polymerase	100 reactions	12361010
	500 reactions	12361050
	2,500 reactions	12361250
Platinum SuperFi II PCR Master Mix	100 reactions	12368010
	500 reactions	12368050
	2,500 reactions	12368250
Platinum SuperFi II Green PCR Master Mix	100 reactions	12369010
	500 reactions	12369050
	2,500 reactions	12369250

Long-range PCR: amplification of 20–40 kb sequences




Figure 3. Amplification of long fragments. A 20 kb target from human gDNA and 30–40 kb targets from *E. coli* gDNA were successfully amplified using Platinum SuperFi II DNA Polymerase.

Benchtop stability: for high-throughput applications




Figure 4. Stability of assembled reactions at room temperature. A 0.5 kb fragment was amplified by setting up reactions and leaving them at room temperature for 0 hr and 24 hr before PCR. Results from using Platinum SuperFi II DNA Polymerase (P) and Q5™ Hot Start High-Fidelity DNA Polymerase from NEB (Q) are shown.



Try before you buy. Scan/click to request free sample.

Find more data at thermofisher.com/platinumsuperfi

invitrogen

For Research Use Only. Not for use in diagnostic procedures. © 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. NEB and Q5 are trademarks of New England BioLabs Inc. PrimeSTAR is a trademark of Takara-Clontech Laboratories, Inc. Merck is a trademark of Merck KGaA. KAPA HiFi HotStart is a trademark of Roche Inc. Qiagen and HotStar are trademarks of Qiagen GmbH. **FLY-12680817 1225**

* No purchase necessary. This promotion is available only to life science professionals 21 years or older in the US (excluding Puerto Rico) and Canada who plan to use the sample within 4 weeks of receipt. Sample requests must be received by December 31, 2026, or before sample supplies are depleted, whichever comes first. Limit one free sample per customer. Please allow 6 to 8 weeks for delivery. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply.